

Functional characterization of Nanodisc-embedded GPCR

Key Fidabio benefits

- Detailed in-solution characterization of receptor-ligand interaction
 - Binding affinity (Kd)
 - Nanodisc sizing and quality check
- Low sample consumption
- Built-in quality control - oligomerization and aggregation check

VERSION 02

Riccardo Marabini, Senior R&D and application Scientists at Fida Biosystems

Introduction

G protein-coupled receptors (GPCRs) are the largest family of membrane proteins and they are involved in many vital physiological cell mechanisms, e.g. hormonal signaling, cell cycle regulation, neurotransmission signaling and odor receptors. GPCRs-targeting drugs represent 27% of the global commercial drugs. Additionally, GPCRs are particularly challenging to work with, given their instability and small-molecule ligands that are difficult to label. Therefore, there is an increasing demand for new technologies to advance discoveries and ultimately new treatments for GPCR related disorders.

This work demonstrates that Flow- Induced Dispersion Analysis (FIDA) can be used to fully characterize the binding between Nanodisc-embedded β 2-adrenergic receptor and Nanobody (Nb80). Further, FIDA has built-in quality control Making it possible to assess the receptor's oligomeric state and the integrity of the sample. FIDA is a new capillary-based technology for measuring binding affinity through absolute and accurate size determinations of analytes in a pressure-driven flow system. The change in apparent size forms the basis for the accurate analysis/determination of binding affinity and protein stability.




Figure 1. Schematic representation of the system analyzed in this work, composed of β 2-adrenergic receptor (in green) embedded in a Nanodisc, and Nanobody Nb80 binding into the G-protein binding pocket.

Material & Methods

FIDA 1 instrument with 480 nm LED fluorescence detection for binding experiments (Fidabio ApS). FIDA standard capillary (i.d.: 75 μm , LT: 100 cm, Leff: 84 cm). HEPES buffer 10mM pH 7.4, 0.5 mg/ml BSA was used as the working buffer. Nanobody Nb80 was used as the indicator. Nb80 was labeled with an Alexa Fluor[®] 488 Protein Labeling Kit from ThermoFisher Scientific.

NanoDisc-embedded β 2-adrenergic receptor (β 2AR) was used as the analyte (0-1 μM). Sample analysis was performed by filling the capillary with the analyte (4 μL), followed by an injection of 39 nL of preincubated indicator+analyte, which was mobilized towards the detector at 400 mbar.

Results

Nb80 binds to the active GPCR

The FIDA technology provides an absolute measurement of the hydrodynamic radius (Rh), and in this work, it was used to measure the size change of Alexa488-labeled Nb80 upon binding to increasing concentration of Nanodiscs-embedded β 2AR (ND- β 2AR). The change in apparent Rh of Nb80 was plotted as a function of increasing ND- β 2AR concentration (0-1 μM) in presence of isoproterenol 10 μM , an agonist that binds to the extracellular side of the GPCR and turns it into the active state (Figure 2A, green line). The Rh of Nb80 without ND- β 2AR was 1,76 nm, and upon addition of ND- β 2AR the size increased up to 5.25 nm, clearly indicating a binding event.

The Rh of 5.25 nm is in agreement with the average size of these Nanodiscs, which is around 5nm. To demonstrate the specificity of the binding, 2 additional controls were measured. The first one was in the presence of an antagonist, alprenolol (APNL) 10 μM , which turns the receptor into the inactive state (Figure 2A, orange dot). The second control was a measure of labeled nanobody in the presence of ND- β 2AR but without agonist nor antagonist (Figure 2A, blue dot).

The two controls Rh at 3.6 and 3.4 strongly demonstrate the reduced binding ability of the GPCR when it is in the inactive state.




Figure 2. (A) Binding curve between Nb80 and ND- β 2AR, IPTR: isoproterenol, APNL: alprenolol. The Rh of Nb80 was plotted as a function of the concentration of ND- β 2AR. **(B)** Raw data showing the indicator peak getting wider as the concentration of ND- β 2AR increased, due to the increase in apparent size upon binding.

Conclusion

The data shows how FIDA can characterize Nanodisc-embedded proteins, in solution and non-invasively. Besides receptor binding, FIDA technology also provides the possibility of assessing

dimerization and oligomerization of GPCRs, and it can work with any scaffold up to 500 nm radius, such as proteoliposomes.